pedalo®

Original Pedalo®

551001, 552001, 552101, 552401,
553011, 553600, 558001, 1853225

www.pedalo.de

Liebe Kundin, lieber Kunde,

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Dear customers,

movement is our element! Our employees occupy since nearly 50 years with the aim to extend to you high quality coordination devices for the daily use in play, sports and therapy. In many loveless produced mass products, the properties functionality and longevity are replaced by "fast broken" and "not repairable". Throw away and buy something new – a trend of this time which we don't expect of our customers.

Quality is a tradition!

We still have them – the solid products.

- Designed with fun and joy and the knowledge of how movement works and for what it helps.

- Made in established tradition in the middle of biosphere territory Swabian Alb.
- Using long lasting and environment-friendly raw materials, we are producing products which "grow again".

It was our aim in the past and always will be it in the future to offer our loyal customers only the best. Therefore, your suggestions are always very welcome and are seen as our challenges for a quality which is maintained in the future.

Ihr Martin Moser (Managing director Holz-Hoerz GmbH)

In 1963

three wheels and two narrow foot boards were the starting signal for a unique company and success story, as it turned out later. We are talking about the unique Pedalo-Sport, which should be used originally as a spontaneous development for therapeutic purposes. Apart from a meaningful therapy device, the Pedalo quickly changed into a play, sports and fun device for everyone.

In 1970

Erich Hoerz, developer of the Pedalo, noticed that many humans were no more able to ride the Pedalo right away. Their balance seems to be getting worse. Therefore, he developed the Pedalo-Classic with longer, not turning platforms and six wheels, which is much easier on the level of difficulty.

In 1985

the Pedalo-Reha ‚S‘ was the next following Pedalo version. Developed especially for the therapy application, it helps those with balance disorders and in rehabilitation after injuries or strokes to recover and improve their capacity.

Made in Germany

Since 1963 from the start of the first brand Pedalo, the company Holz-Hoerz stands for quality products which are 100 % made in Germany. This assures the consumer to get safe products from controlled manufacturing in observance of the highest European safety standards. At the same time we avoid long transport routes to protect environment.

- Made to ecological principles with wood of reforestation
- Proven to be free of prohibited plasticizers
- Tested safety
- Sophisticated in specification and application with 50 years of know-how
- Very stable for children and adults
- 5 years guaranty and lifetime replacement guarantee on the Original-Pedalo
- Difficulty degree alterable for play, sports and therapy
- High gaming value and benefit
- Combinable with other Pedalo-products

Quality warranty

We have high requirements for our products. All products have to motivate young and old as well as people with and without disabilities. All products need to have a high load capacity and a long life expectancy. Our Pedalo-products bear the GS-seal for proven safety. Pedalo-Sport, Pedalo-Classic and Pedalo-Combi have been tested by TÜV Austria and fully comply in load and pollutant test.

Awarded the quality seal 'Bewegte Innovation'

The 'Bewegte Innovation' quality seal of the Würzburg-based movement laboratory InBuB is affixed to appropriately tested and positively assessed products. It is used to determine the movement, learning and play value of movement equipment. Original Pedalo are manufactured in accordance with the Toy Directive 2009/48EG and bear the CE mark. They are tested for stability (EN 71-1) and harmful substances (EN 71-3) in accordance with the standard for toys. They can therefore be used without hesitation.

Pedalo in play, sports & therapy

Already for children from three years, Pedalo are a great challenge full of fun and joy. Besides, they develop their motor skills. Pedalo-products are wonderful to combine with many play and exercise forms. Integrated in ball or partner games, they offer many activities, both indoors and out. Therefore, they contribute to a diversified movement program. With Pedalo, movement is in children's mind and they are really close to it.

Motor objectives:

- Learning of free and auto-nomous Pedalo riding
- Improvement of balance skills
- Expansion of coordinative capacities and development of new movement skills
- Riding the Pedalo for a 10 m distance without any support

Social learning objectives:

- Working cooperatively, offer help one another and also accept the offered assistance
- Learning of consideration
- Rating the own ability and speed
- Applying knowledge in dealing with other persons

Cognitive and affective learning objectives:

- Learning and accepting rules in dealing with others
- Extending knowledge and encourage ambition
- Enhancing safety and self-confidence

Pedalo in play, sports & therapy

Pedalo-products are used for specific exercise sequences in leisure and competitive sports such as soccer, handball, ski jumping, fitness training, golf, volleyball to improve the coordinative skills and to optimize the quickness of action/reaction as well as the performance. Athlete's ability to perform in sports is always the sum of his own coordinative skills which he is able to carry out exactly when needed.

The implementation correlates with the energetic resources and the capacity of the cardiovascular system. All factors of performance ability may be improved through training.

The Pedalo optimizes the total body awareness:

- Interactive combination of movement sequences (capability to couple)
- Response on complex situations and movement processes (responsiveness/ anticipation/ability to adapt)
- Response within spatial conditions (spatial orientation skills)
- Interaction of the postural musculature and motor function (ability to balance = vestibular motor function)

Moreover, it strengthens the upper and lower leg musculature. With an improved vestibular motor coordination, the performance in nearly all kind of sports is increased. The more the training with Pedalo includes different sense stimulations onto the balance, the more safe and efficient the practitioner is going to move.

Since many years the Pedalo is used from therapists in their own practice as well as in rehabilitation and physiotherapy clinics. Through different de-grees in difficulty, Pedalo offers an ideal possibility to adapt on user's current performance capacity. No matter if adults, senior citizens or patients, Pedalo helps many people to maintain and regain their health, well-being and independence.

Pedalo is multi-usable in physiotherapy. It is ideally used in:

- Prophylaxis therapy to prevent injuries in knee or hip area and back pain
- In post-operative therapy of e. g. knee- and total hip endo-protheses (TePs), femoral neck fractures and disk surgeries
- In neurological rehabilitation such as after stroke, Par-kinson, MS
- In pediatrics as training device in a playful way

In physiotherapy it is particularly used in standing position, whereby the trunk erecting and stabilizing muscles are optimally strengthened. To ride Pedalo you need to carry out an alternating steady movement with which the thigh and pelvic muscles as well as the calves are strengthened.

This constant interplay of different muscle groups trains the interactive muscle coordination. Thus, the balance is trained and a dynamic walking movement, respectively an increased feeling of safety when walking or moving in general is imparted. The requested movement sequences when riding Pedalo, foster not only the muscle stability, but also the concentration. Therefore, the cognitive development. Training with Pedalo offers a new movement feeling to the patients and let them feel fun and joy when moving.

The classics

Pedalo are the classics in kindergarten, school, clubs, private use, in sports and therapy. In a playful way they improve coordination, balance and responsiveness. Best quality designed for a daily use with children and adults. Ideally used on even surface as possible, either indoors or outdoors.

EN

The specials

With our special line we meet the manifold requirements of our Pedalo-users. Whether if slalom or riding with seesawing foot plates to get more mobility in the ankle joint.

The Aqua

Pedalo-SAqua are mainly meant for therapeutic applications in water in order to regain the mobility after injuries in locomotor system. Used for play in public pools (in assigned areas) Pedalo-Aqua brings fun and variety.

The Reha-Bar


Featured with side bars, these Pedalo is important aid to regain lost motor abilities after injuries in postural system but also after neurological deficiencies. Due to the individually screwed wheels the Pedalo-Reha-Bar S allows an optimal adjusting on patient's individual needs by additional mounted spacers.

Pedalo® System S

Pedalo S series is characterized by wheels bolted to the axle. The axle provides increased reserves of stability. Designed for use in rehab and therapy for:

severe balance disorders, weak leg muscles, mild hemiplegia, mild hemiplegia after stroke, after surgery especially in the knee and hip area, for spastics even without standing function, wheelchair users with residual strength.

Pedalo® Sport

In leisure and competitive sports the Pedalo-Sport is an important part of training for the optimization of balance, reflexes, posture and movement coordination. Using Pedalo-Sport improves the ankle joint stability and motor ability.

- For ages 5 and up
- Freely rotating foot plates of 6x14 cm each
- Load capacity 120 kg
- 37.5x22x22 cm
- Weight 2.8 kg
- Tires black (551001)
- Accessory: Pedalo-Support Rope (556001-01)

Expert opinion

The Pedalo-Sport is a well-known psychomotor exercise equipment. Where ever you see it, it generates immediately interest, curiosity and the needs to try it right now spontaneously. All age groups likewise enjoy this challenge. Children do not learn riding a Pedalo-Sport at short, it should be continually provided in sports lessons or breaks at schools etc.

Often when the Pedalo-Sport is at free disposal, a couple of children are showing an enormous perseverance and are able to ride only after a few attempts. Whereas others, after initial failures, only achieve notable successes after required support by a partner or with the aid of a rod for example.

My experiences have shown that children prefer to ride forwards, but riding backwards presents a further challenge which should not be neglected particularly.

Alexander Butte, sports teacher

EN

Pedalo® Classic

The Pedalo-Classic is the ideal entry-level model for children and adults. The Pedalo-Classic is equipment that trains the entire body to improve posture and movement coordination. Riding on a Pedalo simulates walking and requires the fluent load change between left and right leg. As a result, the responsiveness, the leg axis stability and thus the safe walking are improved. Support bars can not be attached to this model.

- For age 3 and up
- Foot plates of 30x14 cm each
- Load capacity approx. 200 kg
- 44x37x22 cm
- Weight 6.2 kg
- Tires black (552001)
- Accessory: Pedalo-Foot Pads (144056)
- Conversion into 2 pcs Pedalo-Sport possible

Expert opinion

The Pedalo-Classic fosters the coordinative and perceptual skills and is excellently suited to promote mobility and physical activity. Riding Pedalo with helping, coordinating and arranging among each other and especially playing together, all these have positive effect on social behaviour.

With a Pedalo-Classic you can ride nearly everywhere: At school you ride during sports education in gyms or in the schoolyard during a break. In the leisure time you do at home in the house or outdoors on the footpath. Even there are slight steep and slopes, with a little practice this will be mastered as well.

Alexander Butte, sports teacher

Pedalo® Combi

This Pedalo provides space for 2 people. Therefore, it allows tandem exercises. Partner exercises have a higher difficulty level, need closer coordination and foster the team spirit. To execute the tasks, extreme concentration is required.

- For ages 3 and up
- Foot plates of 60x14 cm each
- Load capacity approx. 200 kg
- 64x37.5x22 cm
- Weight 7.7 kg
- Tires black (552101)
- Accessories: Pedalo-Support Bar (20090000) and Pedalo-Foot Pads (144056)

By using the Pedalo-Combi the partners are faced with new challenges and need quickly adapt to changing circumstances.

Such as riding back to back, face to face or both in the same direction.

Pedalo® Family

Fun, excitement and variety. A pause game device for play-grounds and kindergarten, which motivates groups of up to four people to move. Assimilation, coordination, alignment and integration to each other and with each other are particularly demanded in this case. The two support bars are included for additional stability and safety.

- For age 3 and up
- Foot plates of 100x14 cm each
- Incl. Support Bars Fixed (20090000)
- Load capacity approx. 200 kg
- 104x37.5x22 cm (w/o support bars)
- Weight 13.8 kg
- Tires black (553011)

Riding as a pair or in groups of three or four, requires high attentiveness and consideration of the others. Each individual activity is followed by reaction of the others. The better the team, the better the riding rhythm. Particularly to start and stop riding, a good communication is required. The individual's social and communicative competencies and team-work ability is shown soon.

Pedalo® Reha-Bar S

The special features of this Pedalo are specifically designed for requirements of therapy and senior citizens' sport. It is very stable. Through added spacers between the wheels, as well as through the long adjustable telescopic support bars (including two side bars), it provides a broader and more secure standing area.

- Foot plates of 100x14 cm
- Load capacity approx. 280 kg
- 104x45x75-115 cm (with support bars)
- Weight 21.5 kg
- Tires grey (18553225)
- Accessory: Pedalo-Foot Pads (144056)

Pedalo® Classic S Aqua

The Pedalo-S Aqua is used in hydro therapy to treat balance dysfunctions and responsiveness, to stabilize and mobilize the loco motor system especially that of the legs and the lower back in patients with arthrosis in knee or hip area. With athletes after injuries. With the help of Pedalo-S Aqua a reasonable hydrotherapy with simultaneous relieving of musculoskeletal system can be done.

- Anti-slip foot plates of 30x14 cm each
- PVC foot plates, stainless steel axles
- 43x37.5x22 cm
- Weight 8.5 kg
- Tires blue (552401)

Expert opinion

Aqua pedalo riding in water is more than just a leisure activity that is a lot of fun. Riding in the water offers the advantage, especially for seniors and rehabilitation patients who are at risk of falling, that the movement can be performed while relieving the joints and thus gently, fall-free and without danger. Therefore, the Aqua-Pedalo is often used in therapy after ankle injuries, cruciate ligament injuries, knee and hip operations, intervertebral disc problems and for rheumatism patients to improve balance, coordination and mobility.

The postural and locomotor muscles are coordinated and adapted to the rhythmic high-low movement of walking through the constant interplay between tension and relaxation.

The Aqua-Pedalo adds variety to fall prevention and rehabilitation training.

Annette Jacoby (sports teacher, aqua trainer, lecturer for health sports in the movement space water)

EN

Pedalo® Wawago

This Pedalo-Wawago is a special form of Pedalo and gives you an entirely new feeling of motion. Depending on different standing positions, the toes and heels move up and down just like they do during the natural process of walking. Just turn over one of the boards to lengthen it. It mobilizes the foot joints, resulting in stretching and contracting of the calf and shin muscles. Exercising in squat intensifies this effect. Turn over one of the boards to lengthen it. You will experience a new kind of motion (similar to telemark landing in ski-jumping). Try with right foot in front, afterwards with the left foot in front. In rehabilitation and therapy fields the Pedalo-Wawago helps to prevent phlebitis and effects strengthening of the back muscles through asymmetric stimuli).

- For age 6 and up
- Foot plates of 36x14 cm
- Load capacity approx. 150 kg
- 52x42x22 cm
- Weight 10.0 kg
- Tires black (558001)

Pedalo® Slalom

This is the only Pedalo which allows you to take curves. Dou-bled additional wheels, which are freely rotating in principle, may be tightened individually. Hence, the steering capability can be adjusted from riding only in straight direction up to taking curves in different parameters what requires different steering techniques. Taking curves will be initiated by swinging your upper body or by a swinging movement of your trunk. It is an ideal training device for recreational and professional athletes, much appreciated by skiers and snowboarders. For-wards - backwards - slalom - circle on a point - clockwise -counter-clockwise.

- For age 7 and up
- Rotating foot plates of 6x14 cm
- Load capacity approx. 100 kg
- 37.5x22x22 cm
- Weight 3.6 kg
- Tires black (553600)
- Accessory: Pedalo-Support Rope (556001-01)

For special events

Special events need special products:

By request we are pleased to provide your special Pedalo for your event or your individual application. Get in touch with us, we are happy to advise you about the possibilities.

The world's longest Pedalo

More than 20 people thunder in harmony over the asphalt. A coordinative master performance of all! When the Pedalo has gath-ered speed, there is no way out for the individual. Only a clear message and a joint action may bring the train to stop.

The first attempts

- Always go on the lower foot board first.
- For the beginning take a fixed point at eye-level and move towards.
- Carry out your body weight evenly and fluently from one leg to the other.
- For beginners only do with sturdy shoes.

Exercises with the Pedalo

Exercise aids

- Riding with partner (partner join hands).
- Riding along a wall.
- Riding by help of attached Pedalo-support-ropes.
- Riding with Pedalo-Support Bars (e. g. with senior citizens or people with disabilities).
- Riding between gymnastic boxes.
- Riding between bars or along the bars.
- Riding with stilts (similar to cross-country skiing sticks).
- Move hand over hand along a climbing rope while riding (in gyms).

Riding Pedalo®

Riding in sitting position mobi-lizes the entire body. Only by steady charging and releasing of the legs and the upper part of the body, the Riding-Pedalo will slowly pick up speed. This kind of riding not only brings you in motion but also plenty of laughter. It is also frequently requested for therapeutic use.

Riding exercises

- Riding forwards and backwards.
- Changing view directions during riding: in front, to the left, to the right, ahead, close your eyes while doing so, but keep on going fully concentrated.
- Approaching a chosen aiming point in maximum speed and with spontaneous breaking.
- Approaching a chosen aiming point as near as possible with closed eyes.
- Moving in a low squat position.
- Riding uphill/downhill.
- Riding over a soft floor mat.
- Both feet are positioned on the Pedalo, the front of the body points to the ceiling.

Riding with small sand bags

- A small bag filled with sand, lenses or beans, is balanced on the head, shoulder or the back of the hand.
- Throw a small sand bag in the air and catch it while riding in a constant rhythm.
- Transport a small sand bag from A to B (relay game).
- Transport a small sand bag from A into a bucket. It is up to the rider to decide whether it is better to throw the bag from distance. Hitting is important! Afterwards return pretty quickly to A to take the next little bag.

Variant: Riding with pins

- Pins are positioned on the foot boards for transport from A to B.

Exercises with Hula Hoop tires

- Circle the tires around the hip.
- Circle hula-hoop tires on the arm while Pedalo riding. Try at first with only one, than one for each arm or several tires on one or both arms.
- Riding forwards and backwards.

Exercises with table tennis (TT) bat

Always ride in rhythm. While doing so:

- Balance a TT ball on the TT bat.
- Bounce a TT ball on the bat.
- Repeat the exercises with one bat in each hand.
- Juggle a foot or hand ball on the TT bat and pass it back.

Exercises with ball

Always ride in rhythm forwards or backwards. While doing so:

- Bounce a ball (left hand, right hand).
- Throw a ball constantly against a wall and catch it. Change distance to the wall.
- Take a throw-in over the head against the wall.
- Throw in a basketball basket, goal, gymnastic box, tire...
- Take a ball around your hips.
- Let a ball drop in a bucket or put/take it in on/from a desk while passing.

Juggling

- Juggle balls or cloths

Exercises with ball and partner

Riding Pedalo and doing with the ball:

- Throw/catch directly to the partner.
- Throw/catch indirectly to the partner (across the floor).
- Throw/catch with one hand.
- Take a throw-in over the head to the partner.
- Head the ball back to the partner.
- Play return pass game by foot (the standing leg needs to be positioned always on the lower foot board).
- Serve and spike (volleyball).

Exercises with additional Scooter Boards

- Use the Pedalo as driving system while kneeling or lying on a scooter board.
- Ride the Pedalo and at the same time pull the partner (who is on the scooter board) by a rope.

Partner exercises on one Pedalo each

- Both partners ride in the same time and hold each other. Side by side, one after another.
- Standing face to face and keep a gymnastic ball between the outstretched arms (do not hold the ball but only clamp it with needed pressure).
- Riding while holding one another on two rods.
- Clamp a ball between two rods and transport it (over a distance).

Further exercises

- Each partner stands with one foot on only one board. Now they try to make the Pedalo moving through alternately pressing.
- Time trial over 10 m
- Relay race
- Change Pedalo. Start up, change on next Pedalo, continue riding. Ideal for team relays.

Riding with a resistance band for more posture stability

The band is around the driver's hip. The band's ends are held from the partner or are fixed on a wall mount.

Ride forwards against the increasing resistance, then backwards by yielding the tension softly:

- Ride backwards against the increasing resistance, then forwards by yielding the tension softly.
- Compensate the resistance coming from the side. From left/right (partner runs alongside). The tension can be varied.

The driver holds the band's ends. The middle of the band is held by the partner or fixed on a wall mount:

- Swing your arms in a counter movement to your feet.
- Ride forwards against the increasing resistance, then backwards by yielding the tension softly.
- Ride backwards against the increasing resistance, then forwards by yielding the tension softly.

EN

Balance exercises

While balancing...

- Keep the foot plates on the same level.
- Move the foot plates up and down minimally (play with the balance).
- Change basic position such way, that the feet are in front alternately (always freeze the position for a few seconds).
- Move rhythmically. First do one turn to the front and then back again. Then two turns forwards and backwards, three turns, four turns and so on.

Further spare parts by request to info@pedalo.de

Original Pedalo® System

1. Foot Board 30x14.5 cm for Pedalo-Classic (511030)
2. Foot Board 60x14.5 cm for Pedalo-Combi (511036)
3. Foot Board 100x14.5 cm for Pedalo-Family, -Reha and -Bar (511038)
4. Outer wheel 22 cm, with serigraph 'Pedalo', without tire, for classics (511001)
5. Inner wheel 22 cm, without serigraph 'Pedalo', without tire, for classics (511002)
6. Outer wheel 22 cm, with serigraph 'Pedalo', tire grey included for Reha-Bar S (520016)
7. Inner wheel 22 cm, without serigraph, tire grey included for Reha-Bar S (520248)
8. Tyre \varnothing 22 mm, PVC, free of prohibited Phthalates, black (513001); grey (513006-01)
9. Foot Plate lower 6x14.5 cm (511028) - as replacement part in case the axis has too much clearance. Screws are retained unchanged
10. Screwable (eccentric) axis galvanized Flange (M10) incl. 157.5x21 mm (512020)
11. Screws for axis. For rubber-tired wheels you need four pieces each axis (512110)
12. Foot Plate upper 6x14.5 cm for Pedalo-Sport (511029) - to convert 1 Pedalo-Classic into 2 Pedalo-Sport you need 4 pieces of foot plates

Dear customer,

the Pedalo-product you have chosen is made of natural wood. Differences in color, structure, grain, small included knots or scars, prove the genuineness of our product. They are natural and therefore features of using original wood and no quality reduction.

Please note:

Warning! The packing is not part of the product and have to be disposed. Keep plastic bags and films out of the reach of babies and young children – RISK OF SUFFOCATION!

Warning! Even if this device is also play-equipment for children, the assembling should be carried out or supervised always by an adult (in case product comes not fully mounted).

The Pedalo is to be used with caution, since it requires great skill, so as to avoid falls or collisions causing injury to the user or third parties. In case of doubt first, use the device only with assistance and/or appropriate protective equipment. When used from children, the direct supervision by an adult is recommended.

Do not use the Pedalo in traffic or on public roads.

The Pedalo is only to be used for the purposes described in this booklet and other possible product information/user manual included with the item.

To go on the Pedalo please take the first step always on the lower platform! It must be ensured that feet do not slip off the platforms.

Only use device on level ground with even and firm surfaces. When carrying the Pedalo never touch or take it between the platforms – RISK OF CRUSHING/ENTRAPMENT!

After longer or frequent use, check that the screws are tight and retighten if necessary. Defective and worn parts should be replaced by original replacement parts immediately.

Oiled or lacquered parts may be wiped with a slightly damp cloth.

Protect item from excessive damp and wet. To be stored only in dry rooms. Ensure any storage in smaller rooms, containers, boxes and the like, only in completely dry condition to avoid wood damages and possible functional impairments resulting of it.

Additional note for Pedalo Aqua Line: Only to be used in swimming pools in specific identified areas with flat water and under adult supervision.

Copyright by Holz-Hoerz GmbH,
printed copies also partially are not allowed.

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Subscribe to the newsletter now and
save 10% on your next purchase.

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